


INSTRUCTION FOR VIEWER CONFIGURING (x86,x64,Windows) FOR WORK WITH PACS SERVER

0. Contact your system administrator if your rights to install and configure applications on the local computer are limited. Install the version **1.9.0¹** of application “Inobitec DICOM Viewer”.

1. Start application.

2. In the main menu activate Network->Services. The dialog is appear



3. Turn on DICOM Listener and enter parameters AETitle and Port. AETitle is an arbitrary identifier (name) of the client application, which must be later registered on the PACS server, so that the server "knows" about the application. Enter, for example, LOCAL_SERVICE_AE.

4. Specify the free local port on the computer where the application is launched. To find an available port on your computer, run the command line and type **netstat -a**

```
C:\WINDOWS\system32\cmd.exe
C:\>netstat -a
Active Connections
Proto Local Address           Foreign Address         State
TCP   support:epmap          support:0               LISTENING
TCP   support:microsoft-ds  support:0               LISTENING
TCP   support:1027           support:0               LISTENING
TCP   support:1026           support:0               LISTENING
TCP   support:netbios-ssn    support:0               LISTENING
UDP   support:microsoft-ds  *:*                    *:*
UDP   support:isakmp         *:*                    *:*
UDP   support:1025           *:*                    *:*
UDP   support:4500           *:*                    *:*
UDP   support:ntp            *:*                    *:*
UDP   support:1900           *:*                    *:*
UDP   support:ntp            *:*                    *:*
UDP   support:netbios-ns     *:*                    *:*
UDP   support:netbios-dgm    *:*                    *:*
UDP   support:1900           *:*                    *:*
C:\>
```

All ports (numbers) specified in the column Local Address are used by you computer. Choose any number from 1 to 65535, not included in the presented set and enter it in the

1 Important information is highlighted in bold

field "Port" of the dialog. After entering all the necessary data, the approximate view will be as follows




5. Click the "Apply" button. If somehow, by mistake, the entered port appears to be busy, then the following dialog will be displayed




If the port is free, no messages are output. In our example port 3000 is free, so we use it.

6. In the main menu activate Network->Server list. The dialog is appear



press the  button for creating new connection with PACS server. All fields except the name of new connection will be empty



In the "Server Name" field enter connection name, for example, REMOTE_PACS_SERVER.
You can not enter a server description and leave it blank.

In the “Server type” drop-down list specify **DICOM**.

In the “Server mode” drop-down list specify **C-MOVE**.


In the “Client name (SCU)” field enter AETitle of you application, entered early – **LOCAL_SERVICE_AE**.

In the “Service name (SCP)” field enter **the AETitle of you PACS server** (the example specifies the name of the DCM4CHEE server, but in a specific case, this name must be found in the appropriate person responsible for the PACS server).


In the “Server host name or IP” field enter corresponding data. Usually it is **IP of host**, where PACS server running (the example specifies 192.168.0.7).

In the “Server port” field enter **DICOM port of PACS server** (the example specifies port 11112).


Turn on auto-download images for generating thumbnails so that each series in the list is represented as a thumbnail of the first image in the series. Enable the "Default Server" option if you want this server to be always active in the server list. Click Apply. Dialogue will take the following form



Click the “Apply” button. To test connection with server click . If you see the following message in the dialog box




then the connection to the server is successfully installed and then click OK.
If the server is not available, following message will shown




In this case try ti check **IP address of host** and **DICOM port of PACS server**. You probably made a mistake while entering or the PACS server is not running on the specified computer.

In the case of the next display




you should check the “Server Name (SCP)” field. It must exactly match the **AETitle of your PACS server**.

7. In the administrative part of your PACS server, you need to add **AETitle**, the **DICOM port of the client application service** (viewer) and the **IP of the computer** where the client application is started. In our case AETitle corresponds to **LOCAL_SERVICE_AE**, port - **3000** (see item 3 and item 4). To find your **local IP address**, use the Win + R shortcut and enter cmd in the window that opens, and click OK. In the open command prompt window, type **ipconfig** and press Enter.




Among the information obtained, find the line IPv4-address - this is your local IP.




Enter this IP address along with AETitle and the port listed above, so that the server "knows" about our DICOM service. On many PACS servers, there is a function to test connection with the client-application. Use it to check the client's visibility by the server.

8. To search for studies on the specified PACS server in the viewer in the main

window on the toolbar, click the  button. Next, select the desired PACS server from

the drop-down list on the top right 

Click the “Search” button. From the selected server, studies should be loaded and appear in the corresponding list



Name	ID	Sex	DOB	Mod	# Img	Region	Date and Time	Import date	Specials
				CT	1110		2003-12-01 12:00:00		Specials^10
FhM2hk		0000	1945-02-01	CT	922		2008-08-21 15:01:22		Vasculaire^A
fuQQFud			1931-07-08	CT	1718		2005-04-11 10:33:28		Thorax^1CT
yT3QmG				CT	1054		2005-11-16 15:25:13		CT1 abdomé
FOYKBU				CT	1754		2006-12-07 13:26:39		Cardiac^Car
F063TE				CT	760		2006-10-19 14:37:39		Cardiaque^A
7Thjq7g				CT	557		2005-02-01 12:00:00		CER/CT